



**PHARMA**

# **CHAPTERWISE NOTES**

**Pharmaceutical Analysis**

**Basics of Pharmaceutical Analysis**

# CLASSIFICATION OF ANALYSIS



## QUALITATIVE ANALYSIS

Identify  
*what* substance  
is present



Using litmus paper or  
phenolphthalein to see if  
a solution is acidic or basic



## QUANTITATIVE ANALYSIS

Measure how much  
of a substance  
is present



Acid-base titration to  
determine concentration  
of HCl or NaOH


## Basics of Titration

*Titration is an analytical technique used to determine the concentration of an unknown concentration of an identified analyte by using a standard solution whose concentration is known.*

## Components of Titration

1. **Analyte**—The solution of **unknown concentration** to be analyzed.
2. **Titrant**—The solution of **known concentration** used to react with the analyte.
3. **Burette**—A long, graduated tube to deliver titrant **accurately**.
4. **Conical Flask** – Holds the analyte during titration.
5. **Indicator**—A Substance that changes **colour** at the endpoint.

6. **Clamp Stand**—Holds the burette steady during the experiment.



## Indicators

- Auxiliary substances that change color at (or near) the endpoint.
- Can be synthetic (**phenolphthalein, methyl orange**) or natural (**red cabbage, hibiscus**)

## End Point

- The point in a titration at which the indicator shows a **colour change**, signifying that the reaction is complete.
- Should be as close as possible to the **equivalence point** (where stoichiometric amounts have reacted)


## Equivalence Point

- **Definition:**

The stage in a titration where the **amount of titrant added** is **chemically equivalent** to the amount of analyte present, according to the reaction's stoichiometry

- **Difference from End Point :**

- **Equivalence Point** = Theoretical (exact stoichiometric completion)
- **End Point** = Practical (when indicator changes color)



## CONCENTRATION TERMS

### 1. Molarity (M)

- **Molarity is defined as the moles of solute per volume of solution in litres**
- Mathematically, it is expressed as
- $M = \text{Moles of solute} / \text{volume of solution in litres}$
- The unit of molarity is **mol L<sup>-1</sup>**.

**Problem 1:** Calculate the molarity of the solution containing 5 g of sodium hydroxide in 450 mL of solution

**Solution:** Given: Mass of NaOH = 5 g, Volume of solution = 450 mL = 0.450 L, Molar mass of NaOH = 40 g/mol

### STEP 1: Calculate moles of NaOH

$$\text{Moles} = \text{Mass} / \text{Molar Mass} = 5 / 40 = 0.125 \text{ mol}$$

$$\text{Moles} = 5 / 40 = 0.125$$

### STEP 2: Calculate molarity

$$M = \frac{\text{Moles of solute}}{\text{Volume of solution in liters}}$$

$$M = \frac{0.125}{0.450} \approx 0.278 M$$

## 2. Molality (m)

- Molality is defined as moles of solute per mass of solvent in kilograms**
- Mathematically, it is expressed as
- $m = \text{moles of solute}/\text{mass of solvent in kg}$

## 3. Normality (N)

- Normality is defined as the equivalent of solute per volume of solution in litres.**
- Mathematically, it is expressed as  
**N = Equivalents of solute/volume of solution in L**
- If 40 grams of sodium hydroxide is dissolved in 1 litre of solution, then the normality is found to be 1, and the solution is often termed a 1 normal solution

## SUMMARY

MOLARITY, MOLALITY, AND NORMALITY			
Symbol	Molarity	Molality	Normality
Symbol	M	m	N
Definition	Moles of solute per liter of solution	Moles of solute per kilogram of solvent	Gram equivalents of solute per liter of solution
Unit	mol/L	mol/kg	eq/L
Formula	$M = \frac{\text{Moles of solute}}{\text{Volume of solution}}$	$m = \frac{\text{Moles of solute}}{\text{Mass of solvent in kg}}$	$N = \frac{\text{Gram equivalents of solute}}{\text{Volume of solution}}$ (or $N = M \times \text{Valency factor}$ )
Temperature Dependence	Yes (volume changes with temperature)	No (mass doesn't change with temperature)	Yes (depends on volume)
Example	1 M NaCl 1 mol NaCl in 1 L solution	1 m NaCl 1 mol NaCl in 1 kg solvent	1 N H <sub>2</sub> SO <sub>4</sub> 1 eq (49 g) H <sub>2</sub> SO <sub>4</sub> in 1 L